

Gas Module SE™

State-of-the-art comprehensive vital signs monitoring *and* multi-gas analysis for any surgical environment.

Maximum Power.
Maximum Performance.
Proven Datascope Quality.

- Unique, straight-forward user interface – a hallmark of the Passport® and Spectrum™ series – easily configurable to meet specific requirements. Ideal for multiple users. Requires no extensive training.
- Bright, easy-to-read, user-selectable real-time Fast O₂, CO₂, N₂O, and Agent waveform and numerics display, status messages, and alarms greatly simplify procedural monitoring – even at a distance.
- Breath-by-breath Fast O₂, CO₂, and N₂O monitoring with automatic 5-agent identification and measurement for real-time display of all critical inspired and expired gases.
- Compact, light-weight design ensures convenient use across clinical environments.
- Optional Remote Color Display monitor available for enhanced visibility.
- Passport 2 and Spectrum provide comprehensive vital signs monitoring capability. Quick action keys offer one touch operation.



Passport 2® with Gas Module SE™

Datascope proudly offers the Gas Module SE, the most powerful performance extension to our Passport 2 and Spectrum monitors – already standards of comprehensive vital signs monitoring excellence worldwide. The Gas Module SE delivers state-of-the-art gas monitoring and analysis capabilities: identification and quantification of inspired and expired Fast O₂, CO₂, N₂O, and 5 anesthetic agents – Halothane, Isoflurane, Sevoflurane, Desflurane, and Enflurane. Providing impressive capabilities and highly reliable performance, this leading-edge breath-by-breath multi-gas analyzer was designed to meet the comprehensive anesthesia monitoring requirements of virtually every hospital – whatever its size, specialty, or patient base.

Datascope®

PATIENT MONITORING

Gas Module SE™

So accurate. So reliable. So simple to use.

Our compact, light-weight monitors are perfectly suited to the surgical setting. The simplified cable management system ensures fast, convenient placement on the anesthesia cart.



Large, bright displays, easily visible from a distance, and the simple user interface make this vital signs/multi-gas analysis system the logical choice.



Spectrum and Gas Module SE together offer a large, bright 12.1 inch color display, up to 8 waveforms with auto-adjustable large numerics, 4 invasive pressures and comprehensive gas monitoring.

Datascope Corp.
Patient Monitoring Division
Worldwide Headquarters
800 MacArthur Boulevard
Mahwah, NJ 07430-0619
Tel: 1.800.288.2121
Tel: 201.995.8000
Fax: 1.800.926.4275
E-Mail Address:
pm_sales@datascope.com

www.datascope.com

Datascope® and Passport 2® are U.S. registered trademarks of Datascope Corp.
Gas Module SE™ and Spectrum™ are U.S. trademarks of Datascope Corp.
Subject to change. © 2005 Datascope Corp.

P/N: 0002-08-6005 Rev E