



Thinks big.

NEW **Bis**
VISTA
IN-TOUCH TECHNOLOGY

BIS VISTA™ Monitor Specifications

Product description:	BIS™ (Bispectral Index™) monitoring system for display of processed data and real-time EEG waveforms
Monitor weight:	4.5 lb (2 kg)
Monitor dimensions:	7.5 in wide x 8 in high x 5 in deep (19 cm x 20.3 cm x 12.7 cm)
Display size:	4 in high x 5.25 in wide (10 cm x 13 cm)
Digital output:	USB ports A and B, RS232 serial port
Power requirements:	100-240 VAC, 50-60 Hz, 0.7 ampere max
Electrical safety:	Conforms to: UL 60601-1, IEC 60601-2-26, CAN/CSA-C22.2 #601.1
Battery backup:	45 minutes at full operation Recharge time: 6 hours
Alarms:	Auditory and visual, user-adjustable limits

BISx™ Specifications

Weight:	10 oz (.284 kg) including integral cable
Dimensions:	3.75 in wide x 2.5 in high (9.5 cm x 6.3 cm)
Cable length:	9 ft (2.7 m) integral BISx cable 4.5 ft (1.4 m) from BISx to sensor connector

Ordering Information

Part number	Description
186-0210 (USA)	BIS VISTA Monitoring System includes:
186-0215 (Outside USA)	<ul style="list-style-type: none">• Monitor• Power cord• Pole clamp• BISx• PIC (Patient Interface Cable) – connects BISx to patient



Introducing BIS VISTA...

The monitor packed with power.

BIS VISTA offers you the same proven consciousness assessment of BIS, enhanced within a power-packed monitor.

New, user-friendly features include:

- Full-featured, stand-alone monitor
- Touch screen navigation
- User-configurable display
- 10-, 15-, or 30-second smoothing rates
- DocuBIS™ integrated record keeping
- Increased case storage
- Alarms with adjustable volumes
- Enhanced processing power
- Extended battery life

NEW **BIS**
VISTA

To learn more about BIS and BIS VISTA, please contact us at 1-888-BIS-INDE(X)
(1-888-247-4633), or visit us on the Web at www.aspectmedical.com
Outside USA: TEL +31.30.662.9140

To learn more about BIS education, visit www.BISeducation.com

Copyright © 2006 Aspect Medical Systems, Inc. Aspect, Bispectral Index, the BIS logo, and BIS are trademarks of Aspect Medical Systems, Inc. and are registered in the USA, EU, and other countries.
BIS VISTA and BISx are trademarks of Aspect Medical Systems, Inc.

080-0584 1.00

